

DIY Open Science Workshop: Hardware Guide

Build Your Own Eyetracker and EEG headset with open-source hardware, low-cost electronics, and 3D-printed parts

[Building an EEG Headset](#)

[Parts](#)

[Description of Parts](#)

[Suggested 3D Printing Settings](#)

[Assembly Tools](#)

[Assembly Instructions](#)

[Remove residual support material & print flaws](#)

[Glue the FRAME together](#)

[Mount the OpenBCI BOARD HOLDER](#)

[Assemble the Comfort Nodes](#)

[Insert OCTANUT pieces \(x21\) into frame](#)

[Identify electrode locations](#)

[Connect wiring to OpenBCI](#)

[Assemble electrode mounts \(x8\)](#)

[Insert electrode units into your Ultracortex](#)

[Connect wires to OpenBCI & connect ear clip electrodes \(Reference & Ground\)](#)

[Fasten the BOARD COVER](#)

[Building an Eyetracker](#)

[Parts](#)

[Tools](#)

[Assembly Instructions](#)


[Prepare Webcams](#)

[De-case Cameras](#)

[Remove the camera lens](#)

Building an EEG Headset

This is a modified version from OpenBCI Mark III Nova guide adopted for this workshop.



The Ultracortex is an open-source, 3D-printable headset intended to work with the [OpenBCI](#) system. The OpenBCI EEG board is capable of recording brain activity (EEG), muscle activity (EMG), and heart activity (ECG).

Parts




Figure 2. 3D printed parts for frame and electrode



Figure 3. Internal electrode parts



Figure 4. Wiring




Figure 5. EEG board



Figure 6. Zip Ties

Description of Parts

Spring 1 (x8) — aka "Soft Spring" — [Century Spring Corp. Stock #: U-4](#) — this spring is not as strong as spring 2. We use it to mount the electrodes. The low "Rate" (1.3 lbs/in) makes it more forgiving and comfortable when holding the spiky electrode against your head.

Spring 2 (x5) — aka "Strong Spring" — [Century Spring Corp. Stock #: S-71315S](#) — this spring has a much stronger "Rate" (14 lb/in). This makes it ideal to act as a standoff at locations where you do not have a spiky electrode.

We use small stainless steel screws and hex nuts to fasten the FRI electrodes (listed above) to the 3D-printed electrode mount components & wiring that connects the electrodes back to the OpenBCI board. We used stripped [Gold Cup Electrodes](#) as the cabling, by removing the gold cup with a wire cutter and looping the exposed wire around the screw between the two tightened hex nuts.

Screws (8x) — [Stainless Steel Pan Head Phillips Machine Screw](#), 2-56 Thread, 3/4" Length

Nuts (16x) — [Stainless Steel Hex Nut](#), 2-56 Thread Size, 3/16" Wide, 1/16" High

Screw for Board Case (4x) — #4 [Drive Screw](#)

[Zip tie](#) (10x)

Gold Cup Electrodes (1 bundle) — We strip apart eight electrodes from the [Gold Cup Electrodes](#).

Dry (spiky) electrodes (6x) — [Disposable / Reusable Dry EEG Electrode](#)

Dry (non-spikey) electrodes (2x) — [Disposable / Reusable Cup Wet/Dry EEG Electrode](#)

Ear Clip reference electrodes (2x) — [TDI-430 Silver-Silver Chloride Ear Clip Electrode](#)

EEG Board (1x) — [OpenBCI 32bit Board](#) for 8 electrode channels

Battery Pack (1x) — ~[500mAh lithium ion rechargeable battery pack](#)

Suggested 3D Printing Settings

FRAME_FRONT & FRAME_BACK

Material: PLA

Supports: YES

Raft: hopefully NO (but if supports aren't sticking, try the raft)

Infill: 20%

Layer Height: 0.2mm

Number of Shells: 3

Speed while extruding: 50-70% (slow it down if possible; these parts are detailed)

MECH_PARTS (OCTANUT / BOLT / SPRING_CASING / ELECTRODE HOLDER)

Material: PLA

Supports: NO

Raft: NO

Infill: 20%

Layer Height: 0.2mm

Number of Shells: 3

Speed while extruding: 50-70% (slow it down if possible; these parts are detailed)

BOARD_MOUNT, BOARD_COVER_STANDARD, & BOARD_COVER_ADVANCED

Material: PLA

Supports: NO

Raft: NO

Infill: 20%

Layer Height: 0.2mm

Number of Shells: 3

Speed while extruding: 50-70% (slow it down if possible; these parts are detailed)

Assembly Tools


1. Super Glue
2. coarse flat & circular files (for removing support artifacts)
3. medium sand paper
4. exacto blade (optional)
5. philips head screw driver
6. wire cutters
7. needle-nose pliers (optional)
8. Snippers (optional)



Assembly Instructions

Remove residual support material & print flaws

Use sand paper, a file, and snippers to clean your FRAME and other 3D-printed parts. The most important part of this process is that you thoroughly clean out the frames nodes where you will place your OCTANUT pieces.




Glue the FRAME together

Carefully, glue the FRAME_FRONT and FRAME_BACK together with your super glue. The best way to do this is to place both halves of the frame on a level surface and carefully bring them together. Be sure to be precise; it's VERY difficult to pull the pieces apart once you've put them together.



Mount the OpenBCI BOARD HOLDER

Use four #4 Drive Screws to mount the BOARD HOLDER to the FRAME. Make sure that the orientation of the BOARD HOLDER matches that of the pictures below:



Assemble the Comfort Nodes

In order to fit the OCTANUT pieces into the frame, you will need the Comfort Nodes in order to securely glue the OCTANUTs in place.

You will need the following pieces:

OCTABOLT
Strong Spring
ELECTRODE HOLDER
Comfort insert




First, put the comfort insert into front of the ELECTRODE HOLDER.



Next, insert the spring into the back of the ELECTRODE HOLDER.
![image](../assets/ultracortex-nova-images/spring_electrodeholder.jpg)




Then, slide the OCTABOLT over the top:




Insert OCTANUT pieces (x21) into frame


Note: at each node, you must line up the cut in the OCTANUT with the indentation on the frame as indicated in the picture below.



For each OCTANUT, insert the OCTABOLT so that the OCTANUT is held into the frame. Now you can glue the OCTANUT to the frame like this:





Your Ultracortex should now look like this:



Embed OpenBCI into the Ultracortex

Connect your ~500 mAh lithium ion rechargeable batter to the back of your 32bit OpenBCI Board. Then fold the battery and its wires neatly behind the board before inserting the board into the BOARD_MOUNT. You can then use one of your #4 Drive Screws to secure the OpenBCI Board to the BOARD_MOUNT. Typically you don't need to use any #4 Drive Screws because the BOARD_COVER locks the OpenBCI Board in place, but in this case it's a good idea because you'll want to see where you're connecting your wires.



Identify electrode locations

Before creating your electrode mounts, it's a good idea to think about where you may want to place the electrodes on the Ultracortex FRAME. The placement of the electrode may affect how



long you make the wire between the electrode and where the OpenBCI is mounted, at the back of the frame.

The Ultracortex node locations are based on the [10-20 system](<https://goo.gl/gSwYw>), which is the internationally accepted standard for electrode placement in the context of EEG.

The images below indicates the default 10-20 electrode locations that the OpenBCI Graphical User Interface expects. This application is great for viewing/recording your EEG and can be found in our [OpenBCI_Processing](<http://>) repo. The blue nodes indicate the 8 default 10-20 locations (channels 1-8) of the 32bit Board. The red nodes indicate the default 10-20 locations of channels 9-16 when using the OpenBCI 16-channel R&D Kit.

For the remainder of this tutorial, the blue nodes on the 10-20 system diagram (channels 1-8 of the OpenBCI default settings) will be used. The channel to 10-20 system correlations are as follows:


- * Channel 1 - Fp1
 - * Channel 2 - Fp2
 - * Channel 3 - C3
 - * Channel 4 - C4
 - * Channel 5 - P7
 - * Channel 6 - P8
 - * Channel 7 - O1
 - * Channel 8 - O2



Connect wiring to OpenBCI

Before shortening any of your wires, connect the female header of each wire to the N pin of channels 1-8, as indicated in the image to the lower right. The N pins are the ones closer to the

OpenBCI Board. By default, the OpenBCI Board references these 8 pins with the SRB 2 pin (the bottom SRB pin).



Measure, cut, and strip wires


For each wire, measure the distance between where it is connected to the OpenBCI Board and the node that it will terminate at on the Ultracortex frame. Give yourself 2-3 inches of slack (extra length), because you're going to strip the end of the wire, and you may want the extra slack to zip-tie/tape the wiring to the frame later on.



As mentioned above, if you're using the default OpenBCI electrode configuration, cut the wires so that they terminate at the following nodes:

- * Channel 1 - Fp1
 - * Channel 2 - Fp2

- * Channel 3 - C3
- * Channel 4 - C4
- * Channel 5 - P7
- * Channel 6 - P8
- * Channel 7 - O1
- * Channel 8 - O2


After you're done measuring, cutting, and stripping your wires, disconnect them from the OpenBCI board. You will need the female header free so you can guide the springs and OCTANUT pieces into place in the following step. But remember which wire goes where!





Assemble electrode mounts (x8)

These are the pieces that comprise one full electrode unit:




First, slide a **Machine Screw (2-56 Thread, 3/4" Length)** through the FRI electrode and guide it into the ELECTRODE HOLDER.




Then, slide the screw and electrode into the ELECTRODE HOLDER





Next, twist a *Hex Nut (2-56 Thread Size, 3/16" Wide, 1/16" High)* onto the screw, securing the electrode to the holder.



Next, twist the exposed metal of your wire so that it is tight and clean.




This part is a little bit tricky! Wrap the exposed metal of your wire once around the remaining thread of the screw. You can use your finger nail to pinch the wire down against the nut, while wrapping.




Then, twist your second Hex Nut onto the screw, securing the wire between the two nuts. You can use wire cutters to firmly screw the nut into place.



Slide a **Suggested Spring 1 (aka "Weak Spring")** and then the OCTABOLT around the wire as shown below.



Snap the full electrode unit together, guiding the semi-spherical protrusions on the outside of the SPRING_CASING into the relief cuts of the OCTABOLT.



Voila!

Insert electrode units into your Ultracortex


Once you've assembled all of your electrode units, twist them into their respective Ultracortex frame nodes. If the parts are tough to twist into place by hand, use the OCTATOOL that comes with your kit as shown in the picture below. Over time the resistance between the OCTANUT and OCTABOLT will diminish, and the pieces will twist into place more easily.




Connect wires to OpenBCI & connect ear clip electrodes (Reference & Ground)

Now that all of your electrode units are connected to your Ultracortex frame, reconnect the electrode wires to the OpenBCI board as detailed above. If you have some handy, you may

want to use zip ties to secure the loose wiring to the Ultracortex frame.




Next, connect two touch-proof adapter cables to SRB2 and bottom BIAS pin (white and black respectively in the image below). You can use the touch-proof adapter that comes with the OpenBCI 32bit Board or OpenBCI 16-channel R&D Kit. Then connect your [ear clip electrodes](<http://fri-fl-shop.com/product/td-430-silver-disc-electrode-ear-clip/>) to these touch-proof adapters. These ear clip electrodes serve as the reference and bias (ground with common-mode noise rejection) for your EEG system.



Fasten the BOARD_COVER

Now clip your board cover into place. If you've soldered any of the header rows onto your OpenBCI board, you'll need to use the BOARD_COVER_ADVANCED, as seen in the picture below. Otherwise, you can use the BOARD_COVER_STANDARD with the OpenBCI logo.



Building an Eyetracker

This is a modified version from Pupil Lab's DIY Eyetracking guide adopted for this workshop.

Parts

- Cheap pair of sunglasses
- Playstation Eye
- 2 IR-LEDs
- Exposed (black) film negative
- Aluminum 9-gauge wire
- Wire Ties
- Alligator Clips
- 200 ohms resistor

Tools

- Solder station, wick, flux (in lieu of alligator clips)
- Small philips screwdriver
- Prying tool to help un-case the webcams
- Box cutter

Assembly Instructions

Prepare Webcams

The first step is to modify the cameras so we can use them for eye-tracking.



De-case Cameras

Take out the rubber holders in the back of the webcam casing.







Remove the camera lens

Unscrew the camera lens from the circuit board.